



# BD™ Cardiovascular Fabrics, Surgical Felts and Tapes



# BD Full Line of Cardiovascular Fabrics

Cardiovascular Fabrics					
Product/Indications	Product Code	Size (cm)	Nominal Thickness (mm)	Size (inches)	Units/Box
<b>Knitted Polyester</b>					
<b>BD Sauvage™ Filamentous Fabric</b> Recommended for patch graft angioplasty such as in carotid endarterectomy procedures and ASD/VSD repairs	007942	2.5 x 10.2	0.61	1 x 4	10
	007944	5.1 x 10.2	0.61	2 x 4	10
	007828	10.2 x 10.2	0.61	4 x 4	5
	007829	15.2 x 15.2	0.61	6 x 6	5
<b>BD DeBakey™ Elastic Knit Fabric</b> Used for patch graft angioplasty and repair of intracardiac defects	007830	10.2 x 10.2	0.57	4 x 4	5
	007831	15.2 x 15.2	0.57	6 x 6	5
<b>Knitted PTFE</b>					
<b>BD™ Edwards Outflow Tract Fabric</b> Lower permeability material used for aortic and pulmonary outflow tract patching and intracardiac defects	007835	15.2 x 15.2	0.91	6 x 6	5

Braided Tapes				
Product/Indications	Product Code	Size (mm x cm)	Size (inches)	Units/Box
<b>BD™ PTFE Braided Tape</b>	007913	2.0 x 61.0	$\frac{1}{16}$ x 24	20
	007914	4.0 x 61.0	$\frac{3}{16}$ x 24	20

For the purpose of tissue approximation, as a temporary or permanent ligature, and to retract vessels during a surgical procedure. BD™ PTFE Braided Tapes are packed 2 per sterile pouch, 10 pouches per box, unless otherwise noted.

Surgical Felts					
Product/Indications	Product Code	Size (cm)	Nominal Thickness (mm)	Size (inches)	Units/Box
<b>BD™ Low Porosity PTFE Felt</b> Low porosity helps to control bleeding	007838	10.2 x 10.2	1.85	4 x 4	5
	007839	15.2 x 15.2	1.85	6 x 6	5
<b>BD™ PTFE Felt</b> Used as a buttress for sutures and friable tissue	007975	1.3 x 10.2	1.65	1/2 x 4	10
	007973	2.5 x 2.5	1.65	1 x 1	10
	007976	2.5 x 15.2	1.65	1 x 6	10
	007977	5.1 x 5.1	1.65	2 x 2	10
	007836	10.2 x 10.2	1.65	4 x 4	5
	007837	15.2 x 15.2	1.65	6 x 6	5
<b>BD™ PTFE Felt (Thick)</b> Extra thickness provides added support to tissue	007958	10.2 x 10.2	2.87	4 x 4	5
	007959	15.2 x 15.2	2.87	6 x 6	5

## BD™ Cardiovascular Fabrics

**Indications for Use:** **1)** The knitted and woven polyester fabrics (BD Savage™ Filamentous fabrics, BD DeBakey™ Double Velour fabrics, BD DeBakey™ Elastic Knit fabrics, BD DeBakey™ Woven fabrics) are indicated for use in cardiovascular surgical procedures requiring patch graft angioplasty such as carotid endarterectomy. These fabrics are also indicated for repair of certain intracardiac anomalies such as septal defects. **2)** BD DeBakey™ Woven Fabrics and Edwards Outflow Tract Fabrics are indicated for aortic and pulmonary outflow tract repair and repair of intracardiac defects. Contraindications: Due to the high permeability of the knitted polyester fabrics, they are contraindicated for use in patients requiring prolonged systemic or high dose heparinization except for use in repair of intracardiac defects.

**Warnings:** **1)** These fabrics must be properly preclotted with nonheparinized blood prior to exposure to full arterial pressure to avoid unnecessary bleeding or blood loss. Preclotting is not required when the fabrics are used in the intracardiac position. **2)** While it is common practice to employ moderate doses of intraoperative heparin, care should be taken not to exceed the manufacturer's recommended dose for such procedures. Excessive amounts of heparin can result in bleeding. **3)** Due to variability in patient response to heparinization, it is essential that the adequacy of anticoagulation during surgery and the precision of neutralization at the conclusion of the procedure be closely monitored. Adherence to a strict protocol, determined by each hospital, can prevent excessive bleeding. **4)** As with any cardiovascular fabric, occasional difficulties with hemostasis may occur. In the event that hemostasis is not easily obtained, the surgeon may wish to consider the following: · Investigate for systemic coagulopathy and treat appropriately. · Compression, as necessary. · Additional

sutures and/or pledgets at the anastomosis, as necessary. · Reversal of heparin with protamine sulfate, as necessary. · Utilization of topical coagulation therapy, e.g., thrombin, as necessary. **5)** As with all woven fabric constructions, care must be taken when trimming the fabric to minimize the potential for fraying at cut edges. Caution is highly recommended for heat sealing on all woven patch edges. If the edges of the fabric are not heat sealed, then sutures must be at least 2 mm from the cut edge. **6)** Do NOT expose PTFE fabrics to temperatures greater than 500°F (260°C). PTFE decomposes at elevated temperatures, producing highly toxic decomposition products.

**Precautions:** **1)** STERILE, unless the package is opened or damaged. Single use only. **2)** These products are sterilized by ethylene oxide. DO NOT RESTERILIZE. **3)** Care should be taken when clamping the fabrics to avoid damage to the fibers and preclot. Only atraumatic, shod-type vascular clamps should be used. **4)** Taper point, noncutting needles are recommended so the fabric fibers will not be cut during suturing. **5)** After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

**Potential Adverse Reactions:** Adverse reactions that may occur with the use of these products or with any cardiovascular implant procedure include: · Anastomatic aneurysms · Implant bleeding · Infection · Perioperative hemorrhage · Tissue erosion

### BD™ Surgical Felts and Tapes

Indications for Use: BD™ PTFE and Polyester Surgical Felts are used in various applications for general, vascular and cardiac surgery. They are most commonly used as a patch, a buttress for sutures, and as a material for replacement of segments of the ventricular myocardium after resection.

BD™ PTFE and Polyester Surgical Tapes are indicated for the purposes of tissue approximation, retraction, and as a temporary or permanent ligature.

**Contraindications:** None known.

**Warnings:** **1)** When used as a permanent ligature, care should be taken with surgical tapes to insure that the knots are secure. **2)** Do NOT expose PTFE products to temperatures greater than 500°F (260°C). PTFE decomposes at elevated temperatures, producing highly toxic decomposition products.

**Precautions:** **1)** STERILE, unless the package is opened or damaged. Single use only. **2)** These products are sterilized by ethylene oxide. DO NOT RESTERILIZE. **3)** Care should be taken when clamping surgical tapes to avoid damage to the braided material. Only atraumatic instruments should be used. **4)** Taper point, noncutting needles are recommended when surgical tapes are used as an implant material so the fibers will not be cut during suturing. **5)** Avoid excessive manipulation of the felt surface prior to use. **6)** To assure that felt fibers are cleanly cut and not pulled apart, it is important to use a sharp cutting device or scissors with sharp blades when cutting felt in to smaller pieces. **7)** After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

**Potential Complications:** Potential complications that may occur with the use of these products or with any vascular implant procedure include: · Infection · Tissue erosion

**Please consult package inserts for more detailed safety information and instructions for use.**

벡톤디킨슨코리아(주)  
서울특별시 강남구 테헤란로 142, 16층

bd.com

